

Examples `dshim`
`dshim('method')`
`dshim('help')`

See also *NMR Spectroscopy User Guide*

Related [method](#) Autoshim method (P)
[newshm](#) Interactively create a shim “method” with options (M)
[shim](#) Submit an Autoshim experiment to acquisition (C)
[stdshim](#) Interactively create a shim “method” (M)

ds1sfrq Bandpass filter offset for downsampling (P)

Description For downsampling, selects a bandpass filter that is not centered about the transmitter frequency. In this way, `ds1sfrq` works much like `lsfrq`. If `ds1sfrq` does not exist in the current experiment, add it by entering `addpar('downsamp')`. The command `addpar('downsamp')` creates the digital filtering and downsampling parameters `downsamp`, `dscoef`, `dsfb`, `ds1sfrq`, and `filtfile`.

Values A number, in Hz. A positive value selects a region upfield from the transmitter frequency; a negative value selects a downfield region.

See also *NMR Spectroscopy User Guide*

Related [addpar](#) Add selected parameters to current experiment (M)
[downsamp](#) Downsampling factor applied after digital filtering (P)
[dscoef](#) Digital filter coefficients for downsampling (P)
[dsfb](#) Digital filter bandwidth for downsampling (P)
[filtfile](#) File of FIR digital filter coefficients (P)
[lsfrq](#) Frequency shift of the `fn` spectrum in Hz (P)
[movedssw](#) Set parameters for digital filtering and downsampling (M)
[pards](#) Create additional parameters used by downsampling (M)

dsn Measure signal-to-noise (C)

Syntax `dsn<(low_field,high_field)>:signal_to_noise,noise`

Description Measures the signal-to-noise ratio of the spectrum by first measuring the intensity of the largest peak in the spectral range defined by `sp` and `wp`, and then measuring the noise in the spectral region defined by the position of the two cursors. The noise value returned from `dsn` is not scaled by `vs`. The interrelations between the signal-to-noise ratio, the noise, and peak intensities can be illustrated by comparing `dsn:$sn,$noise` and `peak:$signal`. In this case, `$sn` is equal to $(\$signal / \$noise) / vs$.

Calculate noise by first doing a drift correction on the noise region.
Noise is defined as:

$$noise = 2x \left(\left(\sum_{i=1}^{np} Y_i^2 \right) / np \right)^{\frac{1}{2}}$$

Y_i^2 values are the square of the drift-corrected amplitude and `np` is the number of points in the noise region.

Arguments `low_field` and `high_field` are the upper and lower frequencies of the noise region to be measured. The default is the position of the two cursors.

`signal_to_noise` is the calculated value of signal-to-noise ratio.

`noise` is the noise value measured within the defined spectral region.

Examples
`dsn:$ston`
`dsn(sp+sp,sp+wp-100)`
`dsn(10000,8000):r1`

See also *User Programming*

Related	<code>dres</code>	Measure linewidth and digital resolution (C)
	<code>peak</code>	Find tallest peak in specified region (C)
	<code>sp</code>	Start of plot (P)
	<code>vs</code>	Vertical scale (P)
	<code>wp</code>	Width of plot (P)

dsnmax Calculate maximum signal-to-noise (M)

Syntax `dsnmax<(noise_region)>`

Description Finds the best signal-to-noise in a specified region.

Arguments `noise_region` is the size, in Hz, of the region. The default is the region between the cursors as defined by the parameter `delta`.

Examples
`dsnmax`
`dsnmax(400)`

See also *User Programming*

Related `delta` Cursor difference in directly detected dimension (P)

dsp Display calculated spectrum (C)

Syntax `dsp<(file<,'nodes'>)>`

Description Using the current table of transitions and intensities, `dsp` recalculates the simulated spectrum (using the current value for the linewidth `slw`) and displays the spectrum. `dsp` can only be used after the `spins` program has been run. If only the linewidth `slw` or vertical scale `svs` have been changed, `dsp` can be used to redisplay the spectrum. If a